

## maxXbond: first regeneration system for DNA binding silica matrices

# New solutions allow multiple reuse of valuable material

*Silica matrices are a key technology for the purification of DNA. Today the rapid isolation of pure DNA samples is essential for a variety of molecular biology protocols in research and commercial applications. Products with silica matrices are of high quality and high value. Their major disadvantage is that they can only be used once because after elution substantial amounts of DNA remain attached to the silica matrix and the binding capacity is reduced. To solve this problem AppliChem GmbH, Darmstadt in cooperation with multiBIND GmbH, Cologne developed the first regeneration system for silica matrices. Two innovative solutions remove all nucleic acids and extraneous material from silica matrices and restore the original binding capacity. The system is commercially available under the name maxXbond offering cost savings of approx. 70% with this unique regeneration technology.*

### Regeneration technology for silica matrices

The unique properties of silica matrices for selective DNA binding (Figure 1a) are the basis for all products related to fast and efficient DNA purification. For over two decades, more efficient and application-oriented systems for DNA and RNA have been developed<sup>1,2</sup>. Glass powder ('glass milk'; 'batch procedure') or silica columns allow quick and efficient purification procedures<sup>3</sup>. Principles and problems associated with silica matrices are summarised for mini columns (Figure 1b). Two major tasks for a successful regeneration system are complete removal of residual DNA and restoration of primary binding capacity. Therefore, for a regeneration technology to be useful, the following prerequisites have to be met:


- Quick and easy handling
- Complete removal of all nucleic acids (both free and trapped)
- No damage to the silica matrix
- Complete regeneration of the DNA binding capacity
- Affordability

The two-component maxXbond system with regeneration buffer 1 (RG1) and 2 (RG2) fulfills these requirements. Rapid, efficient

regeneration of DNA binding columns takes only 6 min. (Figure 2). Quantitative analysis of DNA yields demonstrates that the binding capacity of regenerated columns is the same as for new columns (Figure 3a). Analytical agarose gels<sup>2</sup> and PCR analysis<sup>4</sup> verify that regenerated columns are nucleic acid-free (Figure 3b, c). Additional experimental controls for the quality and purity of DNA from regenerated columns are documented in the detailed product information available from your VWR sales specialist. A single DNA binding column can be reused at least 20 times<sup>5</sup>.

Additional considerations in the development and success of the new maxXbond regeneration system are the innovative characteristics of the new solutions:

- All components of maxXbond are bio-degradable, harmless and non-toxic for humans
- No aggressive acids or bases are used so no damage to materials or equipment is observed even after prolonged incubation
- Catalytic and cooperative properties of the maxXbond components cause a rapid and efficient removal or degradation of biological molecules like membrane fragments, proteins and nucleic acids
- Solutions remain active even in the pH range from 6 to 8



The new maxXbond regeneration system can be applied to all commercially available DNA binding columns that contain silica matrices. Preliminary data indicate that any other DNA binding material like glass powder or minerals can also be regenerated by maxXbond. The new product maxXbond is now available to both academic and industrial scientists who seek to optimise their DNA isolation procedures and save a substantial part of their respective costs. The maxXbond regeneration system and solutions have patents pending.



### New product family for efficient regeneration

The first public presentation of maxXbond at Biotechnica 2005 in Hannover, Germany, generated high interest. The positive feedback of the end users initiated the development of new products optimised for the regeneration cycle. The multiple usage of columns is now supported by a new design of columns and collector tubes.

Increased usage of the regeneration system results in a proportional increase in the demand for standard buffer solutions for DNA purification. This prompted the development of maxXmore: a new universal buffer set for all DNA binding columns with silica matrix. This new five-star buffer system has the following advantages over conventional solutions:

- All five buffers are ready to use
- No ethanol has to be added to the wash buffer
- The maxXmore buffer set can be used for all mini columns with silica matrix
- All solutions can be stored at room temperature

New tests demonstrate that maxXbond also regenerates free silica particles for DNA fragment purification and columns with silica matrix for the purification of PCR products. The latest member of the product family is maxXmore PCR, a buffer set with two optimised solutions for the purification of PCR products with regenerated columns.



### Conclusions

The first regeneration system for silica matrices allows substantial cost savings. The high interest of the scientific community and the positive feedback of the primary users have generated many new ideas for future applications. One can expect that maxXbond may be used with other DNA binding materials and protocols. Presently maxXbond is tested for the use with silica matrices applied to genomic DNA preparations, total RNA isolations or messenger RNA purifications, as well as the regeneration of silica particles, magnetic beads and multi-well plates for automated or high-throughput screening. The optimised version of maxXbond for the complete regeneration of silica-based anion exchange columns will be available soon.

### References

1. Birnboim, H.C. & Doly, J. A rapid alkaline lysis procedure for screening recombinant plasmid DNA. *Nucleic Acids Res.* 7, 1513-1522 (1979).
2. Sambrook, J. et al., eds. *Molecular Cloning: A Laboratory Manual*, 2nd ed. (Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York; 1989).
3. Vogelstein, B. & Gillespie, D. Preparative and analytical purification of DNA from agarose. *Proc. Natl. Acad. Sci. USA* 76, 615-619 (1979).
4. Innis, M. A. et al., eds. *PCR Protocols – A guide to methods and applications* (Academic Press, Inc., San Diego, 1990).
5. Esser, K. et al. Nucleic acid-free matrix: regeneration of DNA binding columns. *Biotechniques* 39, 270-271 (2005).

# AppliChem